

ImmunoPrism Sequencing Index Guide

ImmunoPrism Index#	Index Primer Sequence Read
1	ATCAGC
2	CGATGT
3	TTAGGC
4	TGACCA
5	ACAGTG
6	GCCAAT
7	CAGATC
8	ACTTGA
9	GATCAG
10	TAGCTT
11	GGCTAC
12	CTTGTA
13	AGTCAA
14	AGTTCC
15	ATGTCA
16	CCGTCC
17	GTCCGC
18	GTGAAA
19	GTGGCC
20	GTTTCG
21	CGTACG
22	GAGTGG
23	ACTGAT
24	ATTCCT

Note: For low-level multiplexing for sequencing, it is recommended that the following index primers be used:

Pool of 2 samples: Index 6 and 12 or 17 and 23

Pool of 4 samples: Index 16, 17, 19 and 22

Pool of 3 samples: Index 4, 6, and 12 or 13, 17 and 22

Pool of 6 samples: Index 2, 4, 5, 6, 7 and 12